

For Research Use Only

# BAT2L Polyclonal antibody

Catalog Number: 21455-1-AP



## Basic Information

Catalog Number:	21455-1-AP	GenBank Accession Number:	BC002872	Purification Method:	Antigen affinity purification
Size:	300 µg/ml	GeneID (NCBI):	84726	Recommended Dilutions:	IF/ICC 1:50-1:500
Source:	Rabbit	UNIPROT ID:	Q5JSZ5		
Isotype:	IgG	Full Name:	HLA-B associated transcript 2-like		
Immunogen Catalog Number:	AG14690	Calculated MW:	1535 aa, 166 kDa		

## Applications

Tested Applications:	IF/ICC, ELISA	Positive Controls:	
Species Specificity:	human	IF/ICC : HeLa cells,	

## Background Information

PRRC2B is a Novel m6A reader that is expressed and functional in cerebral endothelial cells.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using BAT2L antibody (21455-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).