

For Research Use Only

# C1orf77 Polyclonal antibody

Catalog Number: 21456-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

21456-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14736

GenBank Accession Number:

BC002733

GeneID (NCBI):

26097

UNIPROT ID:

Q9Y3Y2

Full Name:

chromosome 1 open reading frame 77

Calculated MW:

248 aa, 26 kDa

Observed MW:

26-30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse spleen tissue,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Toshie Sakuma	24618812	Viruses	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 21456-1-AP (C1orf77 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.