

For Research Use Only

C1orf77 Polyclonal antibody

Catalog Number: 21456-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

21456-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14736

GenBank Accession Number:

BC002733

GeneID (NCBI):

26097

UNIPROT ID:

Q9Y3Y2

Full Name:

chromosome 1 open reading frame 77

Calculated MW:

248 aa, 26 kDa

Observed MW:

26-30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse spleen tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Toshie Sakuma	24618812	Viruses	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 21456-1-AP (C1orf77 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.