

For Research Use Only

RHBDD3 Polyclonal antibody

Catalog Number: 21457-1-AP



Basic Information

Catalog Number: 21457-1-AP	GenBank Accession Number: BC002705	Purification Method: Antigen affinity purification
Source: Rabbit	GeneID (NCBI): 25807	Recommended Dilutions: WB: 1:1000-1:6000
Isotype: IgG	UNIPROT ID: Q9Y3P4	
Immunogen Catalog Number: AG14737	Full Name: rhomboid domain containing 3	
	Calculated MW: 386 aa, 40 kDa	
	Observed MW: 38 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HepG2 cells, L02 cells
Species Specificity: human	

Background Information

RHBDD3, also named as C22orf3, controls autoimmunity by suppressing the production of IL-6 by dendritic cells via K27-linked ubiquitination of the regulator NEMO.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21457-1-AP (RHBD3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.