

For Research Use Only

ZNF625 Polyclonal antibody

Catalog Number: 21458-1-AP



Basic Information

Catalog Number:

21458-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14758

GenBank Accession Number:

BC007868

GeneID (NCBI):

90589

UNIPROT ID:

Q96127

Full Name:

zinc finger protein 625

Calculated MW:

306 aa, 35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC: 1:50-1:500

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

IF/ICC : PC-3 cells,

Background Information

ZNF625 is a member of the KRAB-ZNF protein family-one of the largest families of transcriptional regulators in mammals.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed PC-3 cells using ZNF625 antibody (21458-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red).