

For Research Use Only

# PLP2 Polyclonal antibody

Catalog Number: 21482-1-AP

Featured Product

2 Publications



## Basic Information

Catalog Number:

21482-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15740

GenBank Accession Number:

BC109066

GeneID (NCBI):

5355

UNIPROT ID:

Q04941

Full Name:

proteolipid protein 2 (colonic epithelium-enriched)

Calculated MW:

152 aa, 17 kDa

Observed MW:

15-17 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IHC 1:20-1:200

IF 1:50-1:500

## Applications

Tested Applications:

IF, IHC, WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

WB : SMMC-7721 cells, HEK-293 cells, COLO 320 cells

IHC : human colon tissue,

IF : HepG2 cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Xiaojian Li	31535380	J Cell Physiol	WB
Hua Bai	32596309	Biomed Res Int	IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21482-1-AP (PLP2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human colon using 21482-1-AP (PLP2 antibody) at dilution of 1:50 (under 10x lens).



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using 21482-1-AP (PLP2 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).