

For Research Use Only

C19orf44 Polyclonal antibody

Catalog Number: 21503-1-AP



Basic Information

Catalog Number:

21503-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15960

GenBank Accession Number:

BC027869

GeneID (NCBI):

84167

UNIPROT ID:

Q9H6X5

Full Name:

chromosome 19 open reading frame

44

Calculated MW:

618 aa, 66 kDa

Observed MW:

66-71 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: L02 cells, mouse liver tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



L02 cells were subjected to SDS PAGE followed by western blot with 21503-1-AP (C19orf44 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.