

For Research Use Only

IL-1F10 Polyclonal antibody

Catalog Number: 21530-1-AP



Basic Information

Catalog Number:

21530-1-AP

Size:

270 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16154

GenBank Accession Number:

BC103966

GeneID (NCBI):

84639

UNIPROT ID:

Q8WWZ1

Full Name:

interleukin 1 family, member 10 (theta)

Calculated MW:

152 aa, 17 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

IL1F10 is a member of the interleukin 1 cytokine family. This cytokine is thought to participate in a network of interleukin 1 family members to regulate adapted and innate immune responses.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 21530-1-AP (IL1F10 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.