

For Research Use Only

CCDC70 Polyclonal antibody

Catalog Number: 21558-1-AP



Basic Information

Catalog Number:

21558-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15600

GenBank Accession Number:

BC069770

GeneID (NCBI):

83446

UNIPROT ID:

Q6NSX1

Full Name:

coiled-coil domain containing 70

Calculated MW:

233 aa, 29 kDa

Observed MW:

29 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, Transfected HEK-293 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 21558-1-AP (CCDC70 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.