

For Research Use Only

MIRO1 Polyclonal antibody, PBS Only

Catalog Number: 21560-1-PBS



Basic Information

Catalog Number:

21560-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15605

GenBank Accession Number:

BC060781

GeneID (NCBI):

55288

UNIPROT ID:

Q8IXI2

Full Name:

ras homolog gene family, member T1

Calculated MW:

618 aa, 71 kDa

Observed MW:

70-80 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21560-1-AP (MIRO1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21560-1-PBS in a different storage buffer formulation.