

For Research Use Only

# ZNF699 Polyclonal antibody, PBS Only

Catalog Number: 21562-1-PBS



## Basic Information

Catalog Number:

21562-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15990

GenBank Accession Number:

BC109268

GeneID (NCBI):

374879

UNIPROT ID:

Q32M78

Full Name:

zinc finger protein 699

Calculated MW:

642 aa, 74 kDa

Observed MW:

74 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human brain tissue were subjected to SDS PAGE followed by western blot with 21562-1-AP (ZNF699 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21562-1-PBS in a different storage buffer formulation.