

For Research Use Only

CYP24A1 Polyclonal antibody

Catalog Number: 21582-1-AP

Featured Product

13 Publications



Basic Information

Catalog Number:

21582-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16134

GenBank Accession Number:

BC109083

GeneID (NCBI):

1591

UNIPROT ID:

Q07973

Full Name:

cytochrome P450, family 24, subfamily A, polypeptide 1

Calculated MW:

514 aa, 59 kDa

Observed MW:

54-59 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

IHC 1:20-1:200

Applications

Tested Applications:

IHC, IP, WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: A549 cells, HeLa cells

IP: HeLa cells,

IHC: human kidney tissue, human prostate cancer tissue

Background Information

CYP24A1(1,25-dihydroxyvitamin D₃) 24-hydroxylase, mitochondrial) catalyzes the NADPH-dependent 24-hydroxylation of calcidiol (25-hydroxyvitamin D₃) and 1-alpha,25-dihydroxyvitamin D₃. Cholecalciferol is metabolised by the 25-hydroxylases (CYP2R1 and CYP27A1) in the liver, before 25-hydroxycholecalciferol (25(OH)D₃) relocates to the circulation (PMID:20362668). The gene CYP24A1 has three isoforms, for the transit peptide "mitochondrial domain" removed, the MW of CYP24A1 will be 55-59 kDa, 47-51 kDa and 43 kDa. This antibody is specific to CYP24A1.

Notable Publications

Author	Pubmed ID	Journal	Application
Thomas Knöpfel	28887454	Sci Rep	WB
Arezoo Daryadel	31628396	Sci Rep	WB
Xue-Qi Chen	36255661	Curr Med Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




WB result of CYP24A1 antibody (21582-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-CYP24A1 transfected A549 cells.




A549 cells were subjected to SDS PAGE followed by western blot with 21582-1-AP (CYP24A1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human kidney using 21582-1-AP (CYP24A1 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human kidney using 21582-1-AP (CYP24A1 antibody) at dilution of 1:50 (under 40x lens).



IP result of anti-CYP24A1 (IP:21582-1-AP, 4ug; Detection:21582-1-AP 1:300) with HeLa cells lysate 2800ug.