

For Research Use Only

MAP3K4 Polyclonal antibody, PBS Only

Catalog Number: 21610-1-PBS



Basic Information

Catalog Number:

21610-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16254

GenBank Accession Number:

BC146770

GeneID (NCBI):

4216

UNIPROT ID:

Q9Y6R4

Full Name:

mitogen-activated protein kinase
kinase kinase 4

Calculated MW:

1608 aa, 182 kDa

Observed MW:

170-180 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 21610-1-AP (MAP3K4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21610-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using MAP3K4 antibody (21610-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 21610-1-PBS in a different storage buffer formulation.