

For Research Use Only

# ZDHHC15 Polyclonal antibody

Catalog Number: 21627-1-AP **1 Publications**



## Basic Information

|   |   |  |
|---|---|--|
| <b>Catalog Number:</b><br>21627-1-AP        | <b>GenBank Accession Number:</b><br>BC103982              | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>200 µg/ml                   | <b>GeneID (NCBI):</b><br>158866                           | <b>Recommended Dilutions:</b><br>WB 1:2000-1:16000           |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q96MV8                              |  |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>zinc finger, DHHC-type containing 15 |  |
| <b>Immunogen Catalog Number:</b><br>AG16152 | <b>Calculated MW:</b><br>337 aa, 39 kDa                   |  |
|   | <b>Observed MW:</b><br>42 kDa                             |  |

## Applications

|  |   |
|--|---|
| <b>Tested Applications:</b><br>WB, ELISA | <b>Positive Controls:</b><br>WB : HEK-293 cells, HepG2 cells, SMMC-7721 cells, HeLa cells |
| <b>Cited Applications:</b><br>IP         |   |
| <b>Species Specificity:</b><br>human     |   |

## Background Information

## Notable Publications

| Author         | Pubmed ID | Journal             | Application |
|----------------|-----------|---------------------|-------------|
| Shanmeizi Zhao | 37769867  | Mol Cell Endocrinol | IP          |

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21627-1-AP (ZDHHC15 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 21627-1-AP (ZDHHC15 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.