

For Research Use Only

ZDHHC15 Polyclonal antibody

Catalog Number: 21627-1-AP

1 Publications



Basic Information

Catalog Number:

21627-1-AP

Size:

200 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16152

GenBank Accession Number:

BC103982

GeneID (NCBI):

158866

UNIPROT ID:

Q96MV8

Full Name:

zinc finger, DHHC-type containing 15

Calculated MW:

337 aa, 39 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:16000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

IP

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, SMMC-7721 cells, HeLa cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shanmeizi Zhao	37769867	Mol Cell Endocrinol	IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21627-1-AP (ZDHHC15 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 21627-1-AP (ZDHHC15 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.