

For Research Use Only

Cyclin B2 Polyclonal antibody

Catalog Number: 21644-1-AP

Featured Product

27 Publications



Basic Information

Catalog Number:

21644-1-AP

Concentration:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16329

GenBank Accession Number:

BC105086

GeneID (NCBI):

9133

UNIPROT ID:

O95067

Full Name:

cyclin B2

Calculated MW:

398 aa, 45 kDa

Observed MW:

51 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:12000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

IHC 1:20-1:200

Applications

Tested Applications:

WB, IP, IHC, ELISA

Cited Applications:

WB, IHC, IP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : A431 cells, DU 145 cells, rat testis tissue, HeLa cells, K-562 cells

IP : A431 cells,

IHC : human cervical cancer tissue, human testis tissue, human colon cancer tissue

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

Cyclin B2 (CCNB2) is a member of cyclin family proteins, which regulate the activities of cyclin dependent kinases (CDKs) and different cyclins function spatially and temporally in specific phases of the cell cycle. Cyclin B2 serves a key role in progression of G2/M transition. Cyclin B2 has been found to be up-regulated in a variety of human cancers, such as adrenocortical carcinoma, breast carcinoma, colorectal adenocarcinoma, pituitary adenoma and gastric cancer. The aberrant expression of Cyclin B2 deregulates spindle checkpoints in the cell cycle and results in chromosomal instability (CIN), one of the signature phenotypes of most cancers. Moreover, serum circulating Cyclin B2 mRNA expression has been found increased in cancer patients and associated with cancer stage and metastasis status. Cyclin B2 can be detected double bands (PMID:1651338).

Notable Publications

Author	Pubmed ID	Journal	Application
Fang Li	36116238	Ecotoxicol Environ Saf	WB
Xiaoyue Zhang	34729106	J Cancer	WB
Zhipeng Qi	34822792	Chem Biol Interact	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 21644-1-AP (Cyclin B2 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.




WB result of Cyclin B2 antibody (21644-1-AP; 1:8000; incubated at room temperature for 1.5 hours) with sh-Control and sh-Cyclin B2 transfected A431 cells.




Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using 21644-1-AP (Cyclin B2 antibody at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using 21644-1-AP (Cyclin B2 antibody at dilution of 1:200 (under 40x lens).



IP result of anti-Cyclin B2 (IP:21644-1-AP, 4ug; Detection:21644-1-AP 1:1000) with A431 cells lysate 2000ug.



Non-treated HeLa and thymidine or nocodazole treated HeLa cells were subjected to SDS PAGE followed by western blot with 21644-1-AP (Cyclin B2 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with GAPDH antibody as loading control.