

For Research Use Only

KRTAP3-1 Polyclonal antibody, PBS Only

Catalog Number: 21673-1-PBS



Basic Information

Catalog Number: 21673-1-PBS	GenBank Accession Number: BC113077	Purification Method: Antigen affinity purification
Concentration: 1 mg/ml	GeneID (NCBI): 83896	
Source: Rabbit	UNIPROT ID: Q9BYR8	
Isotype: IgG	Full Name: keratin associated protein 3-1	
Immunogen Catalog Number: AG16508	Calculated MW: 98 aa, 11 kDa	
	Observed MW: 10-15 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

KRTAP3-1 is a member of the keratin-associated protein (KAP) family. The KAP proteins form a matrix of keratin intermediate filaments which contribute to the structure of hair fibers.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21673-1-AP (KRTAP3-1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21673-1-PBS in a different storage buffer formulation.