

For Research Use Only

DCLK1 Polyclonal antibody

Catalog Number: 21699-1-AP

Featured Product

11 Publications



Basic Information

Catalog Number:

21699-1-AP

Concentration:

600 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16292

GenBank Accession Number:

BC152456

GeneID (NCBI):

9201

UNIPROT ID:

O15075

Full Name:

doublecortin-like kinase 1

Calculated MW:

729 aa, 81 kDa

Observed MW:

48 kDa, 82 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC 1:200-1:800

IF-P 1:50-1:500

Applications

Tested Applications:

WB, IHC, IF-P, IP, ELISA

Cited Applications:

WB, IHC, IF, IP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : rat brain tissue, human brain tissue, mouse brain tissue

IP : rat brain tissue,

IHC : human pancreas cancer tissue,

IF-P : mouse brain tissue,

Background Information

DCLK1 (Serine/threonine-protein kinase DCLK1) is also named as DCAMKL1, DCDC3A, KIAA0369 and belongs to the CAMK Ser/Thr protein kinase family. It is a microtubule-associated kinase that can undergo autophosphorylation and it also has microtubule-polymerizing activity that is independent of its protein kinase activity (PMID: 11124993). It plays a unique role in mitotic spindle integrity during early neurogenesis in radial glial cell proliferation and their radial process stability. DCLK1 is a unique marker for distinguishing tumor stem cells from intestinal normal stem cells (PMID: 23202126). This protein has 4 isoforms produced by alternative splicing with the molecular weight of 82 kDa, 81 kDa, 47 kDa and 48 kDa.

Notable Publications

Author	Pubmed ID	Journal	Application
Kai Lu	31562534	J Gastroenterol	WB
Mahan Si	36210463	Cell Biosci	WB
Lanqing Wang	35522902	Clin Transl Med	WB,IHC,IF,IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




rat brain tissue were subjected to SDS PAGE followed by western blot with 21699-1-AP (DCLK1 Antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.




IP result of anti-DCLK1 (IP:21699-1-AP, 4ug; Detection:21699-1-AP 1:1000) with rat brain tissue lysate 4000ug.




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 21699-1-AP (DCLK1 antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 21699-1-AP (DCLK1 antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 21699-1-AP (DCLK1 antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using DCLK1 antibody (21699-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).