

For Research Use Only

BMT2 Polyclonal antibody

Catalog Number: 21744-1-AP



Basic Information

Catalog Number:

21744-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16476

GenBank Accession Number:

BC114615

GeneID (NCBI):

154743

UNIPROT ID:

Q1RMZ1

Full Name:

chromosome 7 open reading frame 60

Calculated MW:

405 aa, 46 kDa

Observed MW:

46-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, U-87 MG cells, mouse liver, rat brain

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21744-1-AP (C7orf60 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.