

For Research Use Only

PURG Polyclonal antibody

Catalog Number: 21750-1-AP



Basic Information

Catalog Number: 21750-1-AP	GenBank Accession Number: BC106708	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 29942	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9UJV8	
Isotype: IgG	Full Name: purine-rich element binding protein G	
Immunogen Catalog Number: AG15820	Calculated MW: 347 aa, 40 kDa	
	Observed MW: 36-40 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells,
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21750-1-AP (PURG antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.