

For Research Use Only

# NELL1 Polyclonal antibody

Catalog Number: 21783-1-AP **2 Publications**



## Basic Information

<b>Catalog Number:</b> 21783-1-AP	<b>GenBank Accession Number:</b> BC069674	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 170 µg/ml	<b>GeneID (NCBI):</b> 4745	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q92832	
<b>Isotype:</b> IgG	<b>Full Name:</b> NEL-like 1 (chicken)	
<b>Immunogen Catalog Number:</b> AG16542	<b>Calculated MW:</b> 810 aa, 90 kDa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
IHC

**Species Specificity:**  
human, mouse, rat

**Cited Species:**  
human, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Mengmeng Liu	26456257	J Endod	IHC
Tang Ronghua R	23264108	J Mol Histol	IHC

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data