

For Research Use Only

Intersectin 1 Polyclonal antibody

Catalog Number: 21862-1-AP **1 Publications**



Basic Information

Catalog Number: 21862-1-AP	GenBank Accession Number: BC116186	Purification Method: Antigen affinity purification
Size: 850 µg/ml	GeneID (NCBI): 6453	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q15811	
Isotype: IgG	Full Name: intersectin 1 (SH3 domain protein)	
Immunogen Catalog Number: AG16517	Calculated MW: 1721 aa, 195 kDa	
	Observed MW: 150 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : fetal human brain tissue,
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Rui Wei	33014778	Front Oncol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 21862-1-AP (Intersectin 1 antibody at dilution of 1:2000 incubated at room temperature for 1.5 hours.