

For Research Use Only

ZNF541 Polyclonal antibody

Catalog Number: 21873-1-AP



Basic Information

Catalog Number: 21873-1-AP	GenBank Accession Number: BC101052	Purification Method: Antigen affinity purification
Size: 350 µg/ml	GeneID (NCBI): 84215	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q9H0D2	
Isotype: IgG	Full Name: zinc finger protein 541	
Immunogen Catalog Number: AG16192	Calculated MW: 1365 aa, 148 kDa	
	Observed MW: 148 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells,
Species Specificity: human, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21873-1-AP (ZNF541 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.