

For Research Use Only

ZNF541 Polyclonal antibody

Catalog Number: 21873-1-AP



Basic Information

Catalog Number:

21873-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16192

GenBank Accession Number:

BC101052

GeneID (NCBI):

84215

UNIPROT ID:

Q9H0D2

Full Name:

zinc finger protein 541

Calculated MW:

1365 aa, 148 kDa

Observed MW:

148 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21873-1-AP (ZNF541 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.