

For Research Use Only

# ZNF541 Polyclonal antibody, PBS Only

Catalog Number: 21873-1-PBS



## Basic Information

Catalog Number:

21873-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16192

GenBank Accession Number:

BC101052

GeneID (NCBI):

84215

UNIPROT ID:

Q9H0D2

Full Name:

zinc finger protein 541

Calculated MW:

1365 aa, 148 kDa

Observed MW:

148 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21873-1-AP (ZNF541 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21873-1-PBS in a different storage buffer formulation.