

For Research Use Only

SHF Polyclonal antibody

Catalog Number: 21909-1-AP



Basic Information

Catalog Number:

21909-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15083

GenBank Accession Number:

BC101601

GeneID (NCBI):

90525

UNIPROT ID:

B3KTY1

Full Name:

Src homology 2 domain containing F

Calculated MW:

480 aa, 53 kDa

Observed MW:

46-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse skeletal muscle tissue, rat skeletal muscle tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21909-1-AP (SHF antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.