

For Research Use Only

ZNF763 Polyclonal antibody

Catalog Number: 21922-1-AP



Basic Information

Catalog Number:

21922-1-AP

Size:

433 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16440

GenBank Accession Number:

BC105713

GeneID (NCBI):

284390

UNIPROT ID:

Q0D2J5

Full Name:

zinc finger protein 763

Calculated MW:

394 aa, 46 kDa

Observed MW:

46-47 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21922-1-AP (ZNF763 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.