

For Research Use Only

FOXR1 Polyclonal antibody

Catalog Number: 21942-1-AP



Basic Information

Catalog Number:

21942-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16625

GenBank Accession Number:

BC028191

GeneID (NCBI):

283150

UNIPROT ID:

Q6PIV2

Full Name:

forkhead box R1

Calculated MW:

292 aa, 33 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : fetal human brain tissue, mouse brain tissue

IP : mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



fetal human brain tissue was subjected to SDS PAGE followed by western blot with 21942-1-AP (FOXR1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



IP result of anti-FOXR1 (IP:21942-1-AP, 4ug; Detection:21942-1-AP 1:500) with mouse brain tissue lysate 4000ug.