

For Research Use Only

TSPAN31 Polyclonal antibody, PBS Only

Catalog Number: 21987-1-PBS



Basic Information

Catalog Number: 21987-1-PBS	GenBank Accession Number: BC010377	Purification Method: Antigen affinity purification
Source: Rabbit	GeneID (NCBI): 6302	
Isotype: IgG	UNIPROT ID: Q12999	
Immunogen Catalog Number: AG17138	Full Name: tetraspanin 31	
	Calculated MW: 210 aa, 23 kDa	
	Observed MW: 30 kDa	

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Tetraspanin-31 (TSPAN31) belongs to the tetraspanin (TM4SF) family, also known as the transmembrane 4 (TM4) superfamily, which is a family of transmembrane proteins with 33 mammalian members. These proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth, and motility. Tetraspanin-31 is thought to be involved in growth-related cellular processes. The experimentally determined molecular mass of 30 kDa is larger than the calculated molecular mass of 23 kDa, which could be a result of post-translational modifications.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human skeletal muscle tissue slide using 21987-1-AP (TSPAN31 Antibody) at dilution of 1:200 (under 10x lens). This data was developed using the same antibody clone with 21987-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human skeletal muscle tissue slide using 21987-1-AP (TSPAN31 Antibody) at dilution of 1:200 (under 40x lens). This data was developed using the same antibody clone with 21987-1-PBS in a different storage buffer formulation.



Rat heart tissue was subjected to SDS PAGE followed by western blot with 21987-1-AP (TSPAN31 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21987-1-PBS in a different storage buffer formulation.