

For Research Use Only

DDX56 Polyclonal antibody, PBS Only

Catalog Number: 21989-1-PBS



Basic Information

Catalog Number:

21989-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17272

GenBank Accession Number:

BC001235

GeneID (NCBI):

54606

UNIPROT ID:

Q9NY93

Full Name:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 56

Calculated MW:

62 kDa

Observed MW:

62 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 21989-1-AP (DDX56 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21989-1-PBS in a different storage buffer formulation.