

For Research Use Only

PDGFA Polyclonal antibody

Catalog Number: 22011-1-AP

1 Publications



Basic Information

Catalog Number:

22011-1-AP

Concentration:

290 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16340

GenBank Accession Number:

BC146594

GeneID (NCBI):

5154

UNIPROT ID:

P04085

Full Name:

platelet-derived growth factor alpha polypeptide

Calculated MW:

196 aa, 22 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xiaoyu Wu	39638514	Carbohydr Polym	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

