

For Research Use Only

ERG Polyclonal antibody

Catalog Number: 22052-1-AP



Basic Information

Catalog Number:

22052-1-AP

Size:

310 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17046

GenBank Accession Number:

BC040168

GeneID (NCBI):

2078

UNIPROT ID:

P11308

Full Name:

v-ets erythroblastosis virus E26
oncogene homolog (avian)

Calculated MW:

486 aa, 55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

