

For Research Use Only

ZNF718 Polyclonal antibody, PBS Only

Catalog Number:22074-1-PBS



Basic Information

Catalog Number:

22074-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17216

GenBank Accession Number:

BC104027

GeneID (NCBI):

255403

UNIPROT ID:

Q35XZ3

Full Name:

zinc finger protein 718

Calculated MW:

478 aa, 55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

ZNF718 (Zinc Finger Protein 718) is a member of the Krueppel C2H2-type zinc finger protein family, a large class of transcription factors characterized by tandem C2H2 zinc finger motifs that mediate sequence specific DNA binding.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using ZNF718 antibody (22074-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red). This data was developed using the same antibody clone with 22074-1-PBS in a different storage buffer formulation.