

For Research Use Only

CNPY4 Polyclonal antibody, PBS Only

Catalog Number:22086-1-PBS



Basic Information

Catalog Number:

22086-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17285

GenBank Accession Number:

BC032339

GeneID (NCBI):

245812

UNIPROT ID:

Q8N129

Full Name:

canopy 4 homolog (zebrafish)

Calculated MW:

248 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 22086-1-AP (CNPY4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 22086-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using 22086-1-AP (CNPY4 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 22086-1-PBS in a different storage buffer formulation.