

For Research Use Only

# SIDT2 Polyclonal antibody

Catalog Number: 22155-1-AP



## Basic Information

Catalog Number:

22155-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17775

GenBank Accession Number:

BC114522

GeneID (NCBI):

51092

UNIPROT ID:

Q8NBJ9

Full Name:

SID1 transmembrane family, member

2

Calculated MW:

832 aa, 94 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data