

For Research Use Only

SIDT2 Polyclonal antibody, PBS Only

Catalog Number: 22155-1-PBS



Basic Information

Catalog Number:

22155-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17775

GenBank Accession Number:

BC114522

GeneID (NCBI):

51092

UNIPROT ID:

Q8NBJ9

Full Name:

SID1 transmembrane family, member

2

Calculated MW:

832 aa, 94 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data