

For Research Use Only

IRE1/ERN1 Polyclonal antibody

Catalog Number: 22203-1-AP

6 Publications



Basic Information

Catalog Number:

22203-1-AP

Size:

560 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17504

GenBank Accession Number:

BC130405

GeneID (NCBI):

2081

UNIPROT ID:

O75460

Full Name:

endoplasmic reticulum to nucleus
signaling 1

Calculated MW:

977 aa, 110 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, mouse, bovine

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shuqi Yan	39560653	Toxicol Ind Health	WB
Chan-Hung Chen	39056589	Environ Toxicol	WB
Yexiao Yan	38821134	Free Radic Biol Med	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

