

Cytokeratin 17 Polyclonal antibody

Catalog Number: 22230-1-AP

Basic Information

Catalog Number: 22230-1-AP	GenBank Accession Number: BC000159	Purification Method: Antigen affinity purification
Size: 550 ug/ml	GeneID (NCBI): 3872	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q04695	IHC 1:200-1:800
Isotype: IgG	Full Name: keratin 17	IF/ICC 1:20-1:200
Immunogen Catalog Number: AG17575	Calculated MW: 432 aa, 48 kDa	
	Observed MW: 48 kDa	

Applications

Tested Applications:
WB, IHC, IF/ICC, ELISA

Species Specificity:
human, mouse

Note-IHC: suggested antigen retrieval with **TE buffer pH 9.0; (*)** Alternatively, antigen retrieval may be performed with **citrate buffer pH 6.0**

Positive Controls:

WB: HeLa cells, MCF-7 cells

IHC: human bowen disease, human cervical cancer tissue

IF/ICC: HeLa cells,

Background Information

Keratins are a large family of proteins that form the intermediate filament cytoskeleton of epithelial cells, which are classified into two major sequence types. Type I keratins are a group of acidic intermediate filament proteins, including K9-K23, and the hair keratins Ha1-Ha8. Type II keratins are the basic or neutral counterparts to the acidic type I keratins, including K1-K8, and the hair keratins, Hb1-Hb6. Keratin 17 is a type I cytokeratin. It is found in nail beds, hair follicles, sebaceous glands, and other epidermal appendages. It is used as a marker for trauma.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




HeLa cells were subjected to SDS PAGE followed by western blot with 22230-1-AP (Cytokeratin 17 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.




Immunofluorescent analysis of HeLa cells using 22230-1-AP (Cytokeratin 17 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).




Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using 22230-1-AP (Cytokeratin 17 antibody at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using 22230-1-AP (Cytokeratin 17 antibody at dilution of 1:50 (under 40x lens).



Immunohistochemical analysis of paraffin-embedded human bowen disease slide using 22230-1-AP (Cytokeratin 17 antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human bowen disease slide using 22230-1-AP (Cytokeratin 17 antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).