

For Research Use Only

CALCA/Calcitonin Polyclonal antibody, PBS Only

Catalog Number: 22244-1-PBS



Basic Information

Catalog Number:

22244-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17737

GenBank Accession Number:

BC093753

GeneID (NCBI):

796

UNIPROT ID:

P01258

Full Name:

Calcitonin

Calculated MW:

141 aa, 15 kDa

Observed MW:

14 kDa, 5 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Background Information

Calcitonin is a peptide hormone consisting of 32 amino acids, primarily synthesized by the parafollicular cells (C cells) of the human thyroid gland. It plays a crucial role in regulating calcium levels in the blood by decreasing them. Calcitonin opposes the actions of parathyroid hormone, which increases blood calcium levels. Calcitonin works through a G protein-coupled receptor, known as the calcitonin receptor, which predominantly transmits signals via the cAMP and PLC/IP3 pathways. Its primary physiological effects are on osteoclasts and the tubular epithelium of kidneys.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



TT cells were subjected to SDS PAGE followed by western blot with 22244-1-AP (CGRP antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 22244-1-PBS in a different storage buffer formulation.



IP result of anti-CALCA/CGRP (IP:22244-1-AP, 4ug; Detection:22244-1-AP 1:600) with TT cells lysate 945 ug. This data was developed using the same antibody clone with 22244-1-PBS in a different storage buffer formulation.