

For Research Use Only

GCAP2 Polyclonal antibody

Catalog Number: 22298-1-AP



Basic Information

Catalog Number:

22298-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002098

GeneID (NCBI):

2979

UNIPROT ID:

Q9UMX6

Full Name:

guanylate cyclase activator 1B
(retina)

Calculated MW:

23 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : rat retina tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Rat retina tissue were subjected to SDS PAGE followed by western blot with 22298-1-AP (GCAP2 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.