

For Research Use Only

MUC2 Polyclonal antibody

Catalog Number: 22307-1-AP

7 Publications



Basic Information

Catalog Number:

22307-1-AP

Concentration:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002457

GeneID (NCBI):

4583

UNIPROT ID:

Q02817

Full Name:

mucin 2, oligomeric mucus/gel-forming

Calculated MW:

540 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human

Cited Species:

human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ata Ur Rehman	36056976	AMB Express	IHC
Hidayat Ullah	39723032	Food Sci Nutr	WB,IHC,IF
Guanru Qian	39250508	PLoS Pathog	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

