

For Research Use Only

# TACR3 Polyclonal antibody

Catalog Number: 22334-1-AP



## Basic Information

Catalog Number:

22334-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17808

GenBank Accession Number:

BC122551

GeneID (NCBI):

6870

UNIPROT ID:

P29371

Full Name:

tachykinin receptor 3

Calculated MW:

465 aa, 52 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

## Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

IHC : human brain tissue, mouse brain tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human brain slide using 22334-1-AP (TACR3 Antibody) at dilution of 1:50.



Immunohistochemical analysis of paraffin-embedded human brain slide using 22334-1-AP (TACR3 Antibody) at dilution of 1:50.