

For Research Use Only

GATA3 Polyclonal antibody, PBS Only

Catalog Number: 22343-1-PBS

Featured Product



Basic Information

Catalog Number:

22343-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17833

GenBank Accession Number:

BC006793

GeneID (NCBI):

2625

UNIPROT ID:

P23771

Full Name:

GATA binding protein 3

Calculated MW:

443 aa, 48 kDa

Observed MW:

37 kDa and 48 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Trans-acting T-cell-specific transcription factor, also known as GATA-3, belongs to the zinc-finger transcription factor GATA family which drives the differentiation of T helper (Th) 2 cells. GATA factors are expressed in a tissue-specific manner and GATA3 is present in both hematopoietic (e.g., T cells) and non-hematopoietic tissues, including the kidney, central nervous system, skin and mammary gland (PMID: 25148023). GATA-3 plays a central role in Th2 responses both in vitro and in vivo (PMID: 16467870). The full length GATA3 protein contains 444 aa corresponding to the molecular weights of 48,0 kDa, and the truncated protein was detected at approximately 37 kDa (PMID: 24758297).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




WB result of GATA3 antibody (22343-1-AP; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-GATA3 transfected MCF-7 cells. This data was developed using the same antibody clone with 22343-1-PBS in a different storage buffer formulation.



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 22343-1-AP (GATA3 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 22343-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 22343-1-AP (GATA3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 22343-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using GATA3 antibody (22343-1-AP) at dilution of 1:600 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), Alpha Tubulin antibody (66031-1-Ig, Clone: 1E4C11, red). This data was developed using the same antibody clone with 22343-1-PBS in a different storage buffer formulation.