

For Research Use Only

CCR3 Polyclonal antibody

Catalog Number: 22351-1-AP

2 Publications



Basic Information

Catalog Number:

22351-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17926

GenBank Accession Number:

BC110297

GeneID (NCBI):

1232

UNIPROT ID:

P51677

Full Name:

chemokine (C-C motif) receptor 3

Calculated MW:

355 aa, 41 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

IHC 1:100-1:400

Applications

Tested Applications:

IHC, ELISA

Cited Applications:

IHC

Species Specificity:

human

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human tonsillitis tissue, human lymphoma tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Lingwei Bu	39032473	Int Immunopharmacol	IHC
Gongjie Ye	38038167	Oncol Rep	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 22351-1-AP (CCR3 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 22351-1-AP (CCR3 Antibody) at dilution of 1:200 (under 40x lens).