

For Research Use Only

CCR3 Polyclonal antibody

Catalog Number: 22351-1-AP **1 Publications**



Basic Information

Catalog Number: 22351-1-AP	GenBank Accession Number: BC110297	Purification Method: Antigen Affinity purified
Size: 600 µg/ml	GeneID (NCBI): 1232	Recommended Dilutions: IHC 1:100-1:400
Source: Rabbit	UNIPROT ID: P51677	
Isotype: IgG	Full Name: chemokine (C-C motif) receptor 3	
Immunogen Catalog Number: AG17926	Calculated MW: 355 aa, 41 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human tonsillitis tissue, human lymphoma tissue
Cited Applications: IHC	
Species Specificity: human	
Cited Species: human, mouse	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Gongjie Ye	38038167	Oncol Rep	IHC


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 22351-1-AP (CCR3 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 22351-1-AP (CCR3 Antibody) at dilution of 1:200 (under 40x lens).