

For Research Use Only

VEGFR2 Polyclonal antibody

Catalog Number: 22457-1-AP **5 Publications**



Basic Information

Catalog Number:

22457-1-AP

Concentration:

680 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17934

GenBank Accession Number:

BC131822

GeneID (NCBI):

3791

UNIPROT ID:

P35968

Full Name:

kinase insert domain receptor (a type III receptor tyrosine kinase)

Calculated MW:

1356 aa, 152 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human

Cited Species:

human, mouse, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xian Wang	36049428	Phytomedicine	WB,IF
Aini Pan	38016610	Int J Biol Macromol	WB
Wei-Xiang Sun	37480413	Curr Med Sci	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data