

For Research Use Only

# NKX3-2 Polyclonal antibody

Catalog Number: 22491-1-AP

1 Publications



## Basic Information

Catalog Number:

22491-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17581

GenBank Accession Number:

BC111966

GeneID (NCBI):

579

UNIPROT ID:

P78367

Full Name:

NK3 homeobox 2

Calculated MW:

333 aa, 35 kDa

Observed MW:

36-40 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse skeletal muscle tissue, mouse small intestine tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yingyun Yang	39485850	Sci Adv	WB,IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 22491-1-AP (NKX3-2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.