

For Research Use Only

HEATR6 Polyclonal antibody

Catalog Number: 22560-1-AP



Basic Information

Catalog Number:

22560-1-AP

Concentration:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18437

GenBank Accession Number:

BC119813

GeneID (NCBI):

63897

UNIPROT ID:

Q6AI08

Full Name:

HEAT repeat containing 6

Calculated MW:

1181 aa, 129 kDa

Observed MW:

129 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB: 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MCF-7 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells lysate was subjected to SDS PAGE followed by western blot with 22560-1-AP (HEATR6 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.