

For Research Use Only

# SEMA6B Polyclonal antibody

Catalog Number: 22571-1-AP



## Basic Information

Catalog Number:

22571-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16707

GenBank Accession Number:

BC142617

GeneID (NCBI):

10501

UNIPROT ID:

Q9H3T3

Full Name:

sema domain, transmembrane domain (TM), and cytoplasmic domain, (semaphorin) 6B

Calculated MW:

888 aa, 95 kDa

Purification Method:

Antigen Affinity purified

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

