

For Research Use Only

# ZZZ3 Polyclonal antibody

Catalog Number: 22608-1-AP



## Basic Information

Catalog Number:

22608-1-AP

Size:

750 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18362

GenBank Accession Number:

BC035079

GeneID (NCBI):

26009

UNIPROT ID:

Q8IYH5

Full Name:

zinc finger, ZZ-type containing 3

Calculated MW:

903 aa, 102 kDa

Observed MW:

150-160 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells,

## Background Information

Zinc finger, ZZ-type containing 3 (ZZZ3) is a component of the ATAC complex, a complex with histone acetyltransferase activity on histones H3 and H4.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 22608-1-AP (ZZZ3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.